Name and Address:

Organization name Rambam-Health care campus

Address Immunology & Tissue Typing Laboratory

P.O.Box 9602

Haifa 31096, Israel

**Phone** +972-4-8542373 **Fax** +972-4-8295245

E-Mail r katz@rambam.health.gov.il

Date of issue 02.01.2014 Version No. 1 Page No. 1 of: 8

#### **Schedule of Accreditation**

Ite m		Бсоре Туре	Site	Types of Test / Properties Measured		Materials / Sa	mples Tested	Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks
Gr	Group of products: Health, Medicine - Hematology									ואה - המטולוגיה	משפחת מוצרים: בריאות, רכ
Bio	Biological Testing, Cell Typing by FACS Analysis									דם ב- FACS	בדיקות ביולוגיות, איפיון תא
1	A	A	P	Cell Markers - Diagnosis and follow up of: Acute leukemias (AML, ALL) Chronic leukemias, (CLL, CML) Lymphomas, patient's treatment follow up of leukemia, lymphoma and BMT. Identification of lymphocyte subpopulations for diagnosis and follow up of primary and secondary immunodeficiencies (including AIDS) and autoimmune diseases, CFS	Cell Markers - Diagnosis and follow up of: Acute leukemias (AML, ALL) Chronic leukemias, (CLL, CML) Lymphomas, patient's treatment follow up of leukemia, lymphoma and BMT. Identification of lymphocyte subpopulations for diagnosis and follow up of primary and secondary immunodeficiencies (including AIDS) and autoimmune diseases, CFS	Bone marrow, peripheral blood, cell suspension of solid tissues biopsies, CSF, peritoneal and pleural fluids	מח עצם, דם היקפי, תרחיף תאים של ביופסיות רקמות מוצקות, נוזל עמוד שדרה, נוזל פלאורלי ופריטונאלי.		Becton Dickinson FACS calibur or FACS canto		AML=Acute Myelocytic Leukemia, ALL = Acute Lymphocytic Leukemia, CLL = Chronic Lymphocytic Leukemia, CML = Chronic Myelocytic Leukemia, BMT= Bone Marrow Transplantation. CFS= Chronic Fatigue Synd. CSF= CerebroSpinal Fluid. FACS= Fluorescence Activated Cell Sorter

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.

Date of issue 02.01.2014 Version No. 1

## **Schedule of Accreditation**

Ite m	Scope Type	Site	Types of Test / Properties Measured		Materials / Samples Tested		Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks
Gro	oup of	product	ts: Health, Medicine	e - Immunology	& Tissue Typing			וסיווג רקמות	ואה - אימונולוגיר	משפחת מוצרים: בריאות, רפ
Bio	logical '	Testing,	<b>HLA typing and Cros</b>	ss Matching by Mai	nual Immunoassays			אימונוכימיות ידניות	של HLA בשיטות	בדיקות ביולוגיות, סיווג והצלבה
2	A	P	Organ Transplantation C	HLA Typing for Organ Transplantation (Class I)	Lymphocytes	לימפוציטים	Serology typing (Class I)	Microscope		HLA- Human Leukocyte Antigen Biotest plates, InnoTrain plates.
3	A	P		Cross Match for Organ Transplantation	Serum, Lymphocytes	סרום, לימפוציטים	CDC- XM	Microscope		CDC- Complement Dependent Cytotoxicity XM- Crossmatch
4	A	P	Haematopoietic Stem Cell Transplantation: Recipient and related potential donors p	HLA Typing for Haematopoietic Stem Cell Transplantation: Recipient and related potential donors (Class I):	Lymphocytes	לימפוציטים	Serology typing (Class I)	Microscope		Biotest plates, InnoTrain plates
5	A	P	for Disease Diagnosis, f	HLA typing:_Support for Disease Diagnosis, Disease prognosis	Lymphocytes	לימפוציטים	Serology typing (Class I)	Microscope		Biotest plates, InnoTrain plates
6	A	P	Paternity testing P	HLA typing for Paternity testing (court ordered)	Lymphocytes	לימפוציטים	Serology typing (Class I),	Microscope		Biotest plates, InnoTrain plates
7	A	P	screening and s	HLA antibodies screening and identification	Serum	סרום	Solid phase-beads	Luminex		Antibody screening and identification kits by Immucor (Lifecodes)

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.

Date of issue 02.01.2014 Version No. 1

## **Schedule of Accreditation**

Ite m	Scope Type		Types of Test / Pro	pperties Measured	Materials / Sa.	mples Tested	Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks		
Gro	Group of products: Health, Medicine - Immunology & Tissue Typing							ה וסיווג רקמות	ואה - אימונולוגיי	משפחת מוצרים: בריאות, רפ		
	Biological Testing, HLA Typing and Cross Matching, Sequencing Based Molecular Assays						קות ביולוגיות, סיווג והצלבה של HLA בשיטות מולקולריות מבוססות קביעת רצף					
8	A	P	DNA Extraction	,	Peripheral blood Bone marrow, Cord blood, Cell suspension of solid tissues, CVS, Amniotic fluid, Nails and Saliva		Extraction by magnetic beads, Manual Extraction	Maxwell 16		Maxwell 16 Blood Purification kit by Promega QIA amp DNA Blood mini kit by Qiagen		
9	A	P	Post transplantation monitoring	Post transplantation monitoring	DNA	DNA	STR analysis	Thermocycler Sequencer		STR- Short Tandem Repeats SGM plus Kit by ABI		
10	A	P	Paternity testing (court ordered )	בדיקת אבהות (צו בית משפט)	DNA	DNA	STR analysis	Thermocycler Sequencer		STR- Short Tandem Repeats SGM plus Kit by ABI		

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.

## **Schedule of Accreditation**

Ite m	Scope Type	Site	Types of Test / Properties Measured		Materials / Sa	mples Tested	Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks		
Gro	Group of products: Health, Medicine - Immunology & Tissue Typing						משפחת מוצרים: בריאות, רפואה - אימונולוגיה וסיווג רקמות					
Biol Assa	_	Testing,	HLA Typing and C	ross Matching ,PCR	R-Based Molecular			PCR מולקולריות מבוססות	של HLA בשיטות	בדיקות ביולוגיות, סיווג והצלבה		
11	A	P	HLA typing for Organ Transplantation (Class I+II):		DNA	DNA	DNA typing (Class I + ClassII), low and high resolution SSP + SSOP	Thermocycler , Luminex		SSP = Sequence Specific Primers SSOP = Sequence Specific Oligo Probes 1 lambda SSP kits, Unitray Kit by Lifetechologies, Olerup ssp Kit , SSO Kit by Immucor (Lifecodes)		
12	A	P	Haematopoietic Stem Cell Transplantation(Class I+II): Recipient, donor- related potential donors, unrelated donor	HLA Typing for Haematopoietic Stem Cell Transplantation(Class I+II): Recipient, donor- related potential donors, unrelated donor (donor registry),	DNA	DNA	DNA typing (Class I + ClassII) , low and high resolution. SSP + SSOP	Thermocycler , Luminex		1 lambda SSP kits, Unitray Kit by Lifetechnologies, Olerup ssp Kit, SSO Kit by Immucor (Lifecodes)		
13	A	P	typing: Support for Disease Diagnosis,	Class I and/or II HLA typing: Support for Disease Diagnosis, Drug allergy, Disease prognosis,	DNA	DNA	DNA typing (Class I + ClassII) , low and high resolution SSP + SSOP	Thermocycler , Luminex		1 lambda SSP kits, Unitray Kit by Lifetechnologies, Olerup ssp Kit, SSO Kit by Immucor (Lifecodes)		
14	A	Р	HLA typing for Paternity testing (court ordered )	HLA typing for Paternity testing (court ordered )	DNA	DNA	DNA typing classI+II, low resolution (high resolution if needed), SSP + SSOP	Thermocycler, Luminex		1 lambda SSP kits, Unitray Kit by Lifetechnologies, Olerup ssp Kit, SSO Kit by Immucor (Lifecodes)		

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.

Date of issue 02.01.2014 Version No. 1

Page No. 5 of: 8

#### **Schedule of Accreditation**

Ite m	Scope Site Types of Test / Properties Measured Materials / San			mples Tested	Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks			
Group of products: Health, Medicine - Immunology & Tissue Typing								וסיווג רקמות	ואה - אימונולוגיו	משפחת מוצרים: בריאות, רפ	
Biol	ogical [	Гesting,	HLA typing and cro	ss matching by FAC	S	בדיקות ביולוגיות, סיווג והצלבה של HLA ב- FACS					
15	A		Cross Match for Organ Transplantation		Serum, Lymphocytes	סרום, לימפוציטים		Becton Dickinson FACS calibur Becton Dickinson FACS canto		Flow XM- Flow Cytometry Crossmatch	

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.

## **Schedule of Accreditation**

Ite m	Scope Type	Site	Types of Test / Properties Measured		Materials / Samples Tested		Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks			
Gro	Group of products: Health, Medicine - Immunology & Tissue Typing						משפחת מוצרים: בריאות, רפואה - אימונולוגיה וסיווג רקמות						
Biol	ogical '	Testing,	Autoimmune Disease	e, Manual Immunoas	ssays			מונוכימיות ידניות	ימוניות בשיטות אי	בדיקות ביולוגיות, מחלות אוטוא			
16	A	P	Anti Epithelial Ab	Anti Epithelial Ab	Serum	סרום	IFA	Fluorescent Microscope		IFA= Indirect Immunofluorescence Kit by INOVA			
17	A	P	ANA (Anti –Nuclear Ab)	ANA (Anti –Nuclear Ab)	Serum	סרום	IFA	Fluorescent Microscope		Kit by ALPHADIA			
18	A	P	APCA (Anti parietal Cell Ab)	APCA (Anti parietal Cell Ab)	Serum	סרום	IFA	Fluorescent Microscope		AMA, ASMA, APCA, are done on the same slide Kit by Kallestad			
19	A	P	AMA (Anti Mitochondrial Ab)	AMA (Anti Mitochondrial Ab)	Serum	סרום	IFA	Fluorescent Microscope		AMA,ASMA,APCA, are done on the same slide Kit by Kallestad			
20	A	P	ASMA (Anti Smooth Muscle Ab)	ASMA (Anti Smooth Muscle Ab)	Serum	סרום	IFA	Fluorescent Microscope		AMA, ASMA, APCA, are done on the same slide Kit by Kallestad			
21	A	P	EMA (Anti Endomysial Ab)	EMA (Anti Endomysial Ab)	Serum	סרום	IFA	Fluorescent Microscope		Slide by INOVA			
22	A	P	GBM (Anti Glomerular Basement Membrane Ab)	GBM (Anti Glomerular Basement Membrane Ab)	Serum	סרום	IFA	Fluorescent Microscope		Slide by INOVA			
23	A	P	ANCA (Anti- Neutrophil Cytoplasmic Ab)	ANCA (Anti- Neutrophil Cytoplasmic Ab)	Serum	סרום	IFA	Fluorescent Microscope		Screening for p-ANCA or c-ANCA Kit by INOVA			
24	A	P	MPO (Anti Myeloperoxidase Ab)	MPO (Anti Myeloperoxidase Ab)	Serum	סרום	ELISA	Elisa Reader		Ab=Antibodies ELISA=Enzyme - Linked Immuno Sorbent Assay Kit by INOVA			
25	A	P	PR-3 (Anti Serine Protease 3 Ab)	PR-3 (Anti Serine Protease 3 Ab)	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA			

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.

Date of issue 02.01.2014 Version No. 1

Page No. 7 of: 8

#### **Schedule of Accreditation**

Ite m	Scope Type	Site	Types of Test / Properties Measured		Materials / Samples Tested		Test principle / Characteristics	Equipment used / Work station	Opinion and Interpretation	Remarks
Gro	Group of products: Health, Medicine - Immunology & Tissue Typing							ה וסיווג רקמות	ואה - אימונולוגיו	משפחת מוצרים: בריאות, רפ
Biol	ogical	Testing,	Autoimmune Disease	, Manual Immunoas	ssays			מונוכימיות ידניות	ימוניות בשיטות אי	בדיקות ביולוגיות, מחלות אוטוא
26	A	P	ASCA (Anti Saccharomyces Cerevisiae Ab)	ASCA (Anti Saccharomyces Cerevisiae Ab)	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA
27	A	P	Anti RNP/SM Ab	Anti RNP/SM Ab	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA
28	A	Р	Anti SM Ab	Anti SM Ab	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA
29	A	Р	Anti SS-A Ab	Anti SS-A Ab	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA
30	A	P	Anti SS-B Ab	Anti SS-B Ab	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA
31	A	Р	`	SCL-70 (Anti scleroderma 70)	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA
32	A	P	tTG (Anti Tissue Transglutaminase Ab)	tTG (Anti Tissue Transglutaminase Ab)	Serum	סרום	ELISA	Elisa Reader		Kit by INOVA

Site: P or T or M, P-Permanent, T-Temporary, M-Mobile

Type of Scopes: A- Fixed, C- Flexible scope in analytical tests: Type of matrix, analytes, experimental systems and/or analytical characteristics may be subject to changes, in accordance with the laboratory's approved and documented procedures. For details, please refer to the list of Accredited Tests, available from the laboratory upon request.